

Alexandrite Laser

Dual-wavelength laser technology | Model LSA-B



BEAURAY-BRANDED PLATFORM

Technology designed for professional settings

A dual-wavelength platform combining 755 nm Alexandrite and 1064 nm Nd:YAG technologies. The platform supports configurable professional treatment workflows and Beauray private-label presentation.

INTENDED PROFESSIONAL APPLICATIONS

- Professional hair-removal applications
- Vascular and skin-rejuvenation applications
- Nail-fungus applications where permitted and appropriately regulated

Platform overview

WAVELENGTHS	755 nm / 1064 nm	COOLING	Compressor + water + air
SPOT SIZE	6-22 mm	CONSTRUCTION	Metal
PULSE WIDTH	5-100 ms	POWER SOURCE	3680 W
DISPLAY	10.4-inch TFT	WARRANTY	2 years

COMMERCIAL & DISTRIBUTION INFORMATION

Manufactured by Sumstars for Beauray branding. OEM/ODM, configuration, language, documentation, availability, and market-specific regulatory requirements are reviewed per project. No regulatory clearance or approval is represented unless explicitly documented.