



Bioteque

Empowering Women's Healthcare

About

Bioteque was founded on a vision to make high-quality, cost-effective OB/GYN solutions accessible to healthcare professionals worldwide. From its earliest days, the company built its foundation on manufacturing excellence in Taiwan, where precision engineering, strict quality controls, and consistent performance define its operations, while also establishing efficient logistics and close, long-standing relationships with healthcare providers.

Today, Bioteque's global headquarters and primary manufacturing operations are based in Taiwan, supported by established operations in the United States, Europe, and the Philippines. This international presence strengthens worldwide distribution and ensures responsive, reliable support across the regions served.

Focused on quality, reliability, and long-term value, Bioteque continues to expand its global reach while delivering dependable OB/GYN devices, accessories, and surgical solutions for diverse clinical environments.

Mission Statement

*For more than 30 years, Bioteque's mission has been simple:
—— To Empower Women's Healthcare.*

This guiding purpose shapes every decision, every innovation, and every partnership. Through trust, quality, regulatory excellence, and continuous progress, Bioteque remains dedicated to supporting healthcare professionals and improving the quality of care for women worldwide—today and for generations to come.




TABLE OF CONTENTS

Pessaries	3-6
Fitting Sets	7
Pessary Accessories	8
Endometrial Samplers	9
Vaginal/Cervical Dilators	10
Bioteque Catheters	11
Electrosurgical Devices	12
Electrosurgical Accessories	13-14
Instruments - Biopsy Punches	15
Instruments - Forceps / Scissors	16
Speculum Series	17-19
Stainless Steel Curettes	19
Stainless Steel Dilators	20
Reference Guide and Selection Chart	21
Odering and Returns	22

Pessaries

Bioteque offers a complete range of medical-grade silicone pessaries in a variety of shapes and sizes to accommodate different patient needs. Each product is manufactured under strict quality control systems and complies with FDA, cGMP, CE, QMS, MDD, ARTG, and ISO 13485 standards, ensuring safety, consistency, and reliability.



Ring

The Ring pessary, available with or without a supporting membrane, is primarily used for mild first- degree uterine prolapse and cystocele.

SKU (W/O support)	R#0	R#1	R#2	R#3	R#4	R#5	R#6	R#7	R#8	R#9
Dimensions (mm)	44	50	57	63	69	75	82	88	95	101
Dimensions (inch)	1.75	1.95	2.25	2.50	2.70	2.95	3.25	3.45	3.75	4.00
SKU (with support)	Rs#0	Rs#1	Rs#2	Rs#3	Rs#4	Rs#5	Rs#6	Rs#7	Rs#8	Rs#9



Ring with Knob

The Ring with Knob pessary, available with or without a supporting membrane, is primarily used for mild first-degree uterine prolapse, cystocele, and stress urinary incontinence.

SKU (W/O support)	RK#0	RK#1	RK#2	RK#3	RK#4	RK#5	RK#6	RK#7
Dimensions (mm)	44	50	57	62.5	69	75	82	87.5
Dimensions (inch)	1.75	1.95	2.25	2.45	2.70	2.95	3.25	3.45
SKU (with support)	RKs#0	RKs#1	RKs#2	RKs#3	RKs#4	RKs#5	RKs#6	RKs#7



Incontinence Ring

The Incontinence Ring is primarily used for stress urinary incontinence, mild first-degree uterine prolapse, and cystocele.

SKU	IRO	IR1	IR2	IR3	IR4	IR5	IR6	IR7	IR8	IR9
Dimensions (mm)	44	51	57	64	70	76	83	89	95	102
Dimensions (inch)	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00



*Bioteque Exclusive

Oval

The Oval pessary, available with or without a supporting membrane, is primarily used for mild first-degree uterine prolapse and cystocele. The pessary’s shape helps to limit rotation.

SKU(W/O support)	OV#1	OV#2	OV#3	OV#4	OV#5	OV#6	OV#7	OV#8	OV#9
Dimensions (mm)	50.5	57	63.5	69.5	76	83	89	96	101
Dimensions (inch)	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.80	4.00
SKU (with support)	OVs#1	OVs#2	OVs#3	OVs#4	OVs#5	OVs#6	OVs#7	OVs#8	OVs#9



*Bioteque Exclusive

Marland

The Marland pessary, available with or without a supporting membrane, is primarily used for second or third-degree uterine prolapse, cystocele, rectocele. The pessary's elevated platform may also relieve stress urinary incontinence.

SKU (W/O support)	M#2	M#3	M#4	M#5	M#6	M#7	M#8
Dimensions (mm)	58	64	70	76	82	88	94
Dimensions (inch)	2.30	2.50	2.75	3.00	3.25	3.45	3.70
SKU (with support)	Ms#2	Ms#3	Ms#4	Ms#5	Ms#6	Ms#7	Ms#8



Gehrung

The malleable Gehrung pessary with a supporting membrane is primarily used for significant cystocele or rectocele, and second or third-degree uterine prolapse. The pessary can be manually shaped to accommodate individual anatomies. Use caution and follow doctor's instructions before undergoing MRI.

SKU	GH#0	GH#1	GH#2	GH#3	GH#4	GH#5	GH#6	GH#7	GH#8	GH#9
Dimensions (mm)	50	55	60	65	70	75	80	85	90	95
Dimensions (inch)	1.95	2.15	2.35	2.55	2.75	2.95	3.15	3.35	3.55	3.75



Shaatz

The Shaatz pessary is primarily used for second or third-degree uterine prolapse and cystocele. Its convex shape ensures a snug fit drainage holes aid in removal.

SKU	SH#0	SH#1	SH#2	SH#3	SH#4	SH#5	SH#6	SH#7	SH#8	SH#9
Dimensions (mm)	38	45	51.5	57	63	70	76	82	88	95
Dimensions (inch)	1.50	1.75	2.05	2.25	2.50	2.75	3.00	3.25	3.45	3.75



Donut

The Donut pessary is primarily used for second or third-degree uterine prolapse, cystocele, and rectocele. This pessary, often used by post-menopausal women, is easily compressed for insertion and removal.

SKU	D#0	D#1	D#2	D#3	D#4	D#5	D#6	D#7
Dimensions (mm)	50	59	65	69	76	82	88	95
Dimensions (inch)	1.95	2.30	2.55	2.70	3.00	3.25	3.45	3.75

Bioteque Pessary CPT Code: A4562 Procedure Code: 57160



*Bioteque Exclusive

FlexiShelf

The Flexishelf pessary is primarily used for second or third-degree uterine prolapse. It has a unique thin curved-base design, rigid inner layer and soft outer layer. The curved cone-shaped handle aids in removal.

SKU	FS06	FS08	FS10	FS12
Dimensions (mm)	57	64	70	76
Dimensions (inch)	2.25	2.50	2.75	3.00



Cube

The Cube pessary, available with or without drains, is primarily used for cystocele, rectocele, second or third-degree uterine prolapse, and complete procidentia.

SKU (w/o drains)	CU#0	CU#1	CU#2	CU#3	CU#4	CU#5	CU#6	CU#7	CU#8	CU#9	CU#10
Dimensions (mm)	25	29	33	37	41	45	50	56	63	70	75
Dimensions (inch)	1.00	1.15	1.30	1.45	1.60	1.75	1.95	2.20	2.50	2.75	2.95
SKU (with drains)	CUd#0	CUd#1	CUd#2	CUd#3	CUd#4	CUd#5	CUd#6	CUd#7	CUd#8	CUd#9	CUd#10



Gellhorn

The Gellhorn pessary with drains is primarily used for cystocele, rectocele, and second or third-degree uterine prolapse. This pessary's stem with knob is easily folded to aid with insertion and removal.

SKU (with drains)	Gd#0	Gd#1	Gd#2	Gd#3	Gd#4	Gd#5	Gd#6	Gd#7	Gd#8	Gd#9	Gd#10
Dimensions (mm)	38	45	51.5	57	63	70	76	82	88	95	101.6
Dimensions (inch)	1.50	1.75	2.05	2.25	2.50	2.75	3.00	3.25	3.45	3.75	4.00



Short Stem Gellhorn

The short-stem Gellhorn pessary with drains is primarily used for cystocele, rectocele, and second or third-degree uterine prolapse. It provides a better fit for some women,

SKU (with drains)	GS#0	GS#1	GS#2	GS#3	GS#4	GS#5	GS#6	GS#7	GS#8	GS#9
Dimensions (mm)	38	44	51	57	63.5	70	76	83	89	95
Dimensions (inch)	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75



Dish

The Dish pessary, available with or without a supporting membrane, is primarily used for second or third-degree uterine prolapse, cystocele, rectocele and stress urinary incontinence.

SKU (W/O support)	DSH#0	DSH#1	DSH#2	DSH#3	DSH#4	DSH#5	DSH#6	DSH#7
Dimensions (mm)	48.5	53.5	59	64	69	75	79	84
Dimensions (inch)	1.90	2.10	2.30	2.50	2.70	2.95	3.10	3.30
SKU (with support)	DSHs#0	DSHs#1	DSHs#2	DSHs#3	DSHs#4	DSHs#5	DSHs#6	DSHs#7



Cup

The Cup pessary, available with or without a supporting membrane, is primarily used for first or mild second-degree uterine prolapse. This pessary is similar to the dish without a knob.

SKU (W/O support)	CP#0	CP#1	CP#2	CP#3	CP#4	CP#5	CP#6	CP#7
Dimensions (mm)	48.5	53.5	59	64	69	74	79	84
Dimensions (inch)	1.90	2.10	2.30	2.50	2.70	2.90	3.10	3.30
SKU (with support)	CPs#0	CPs#1	CPs#2	CPs#3	CPs#4	CPs#5	CPs#6	CPs#7



Hodge

The malleable Hodge pessary, available with or without a supporting membrane, is primarily used for first or mild second-degree prolapse and cystocele. This pessary, is usually prescribed for reproductive-aged women who prefer conservative therapy. Use caution and follow doctor's instructions before undergoing MRI.

SKU (W/O support)	HD#0	HD#1	HD#2	HD#3	HD#4	HD#5	HD#6	HD#7	HD#8	HD#9
Dimensions (mm)	65	70	75	80	85	90	95	100	105	110
Dimensions (inch)	2.55	2.75	2.95	3.15	3.35	3.55	3.75	3.95	4.15	4.35
SKU (with support)	HDs#0	HDs#1	HDs#2	HDs#3	HDs#4	HDs#5	HDs#6	HDs#7	HDs#8	HDs#9

Silicone Inflatable Donut



*Bioteque Exclusive

The Inflatable donut is primarily used for mild second or third degree cystocele in combination with procidentia/prolapse.

SKU	Dimensions	Dimensions
ID 2	51mm	2.00in
ID 3	57mm	2.25in
ID 4	64mm	2.50in
ID 5	70mm	2.75in



Ring Fitting Set
SKU : FS1000



Cube Fitting Set
SKU : FS2000



Dish Fitting Set
SKU : FS3000



Oval Fitting Set
SKU : FS4000



Gellhorn Fitting Set
SKU : FS5000



Ring with support Fitting Set
SKU : FS6000



Ring with Knob Fitting Set
SKU : FS7000



Marland Fitting Set
SKU : FS8000



Cup Fitting Set
SKU : FS9000

Pessary Fitting Set

The Pessary Fitting Set is designed for the physician to determine the proper size and type of pessary for the patient before prescription. 100% Silicone blue in color for identification, Available in 11 types; each set includes 6 sizes, Autoclavable and can be resterilized up to 150 times.

SKU	Type	Size
FS1000	Ring	Size 2-7
FS2000	Cube	Size 2-7
FS3000	Dish	Size 2-7
FS4000	Oval	Size 2-7
FS5000	Gellhorn	Size 2-7

SKU	Type	Size
FS6000	Ring w/ Support	Size 2-7
FS7000	Ringw/ Knob	Size 2-7
FS8000	Marland	Size 2-7
FS9000	Cup	Size 2-7
FS1100	FlexiShelf	Size 6-12



Bioteque Pessary Remover

Designed to aid in removal. It is made from medical-grade stainless steel and can be used with the following types of pessaries:

Ring, Ring with Knob, Oval, Dish, Hodge, Marland and Gellhorn.

Lubricant is required while using the remover. Healthcare professional use only.

SKU	Description	Length	
PRM2	Reusable Pessary Remover	182mm	25units/box

P.O.P.s Ruler

The POP-Q score is a clinical tool used to identify the severity of pelvic organ prolapses and the anatomical success of conservative therapy and surgery.


An easy-to-use measuring device with 0.5cm markings helps perform accurate POP-Q scoring.

SKU	Description	Length	
BA222-222	Disposable P.O.P.s Ruler	180mm	50units/bag

POP-Q QUICK REFERENCE GUIDE

The POP-Q assessment tool measures nine points (see the chart). There are six points measured at the vagina with respect to the hymen. The prolapsed organs are measured in centimeters to the hymen. The measurements are taken when the Valsalva maneuver is performed by the patient lying in the dorsal lithotomy position. Points above the hymen are negative numbers, points below the hymen are positive numbers.

The further three points are the genital hiatus (gh: measured from middle of external urethral meatus to posterior midline hymen), the perineal body (pb: measured from posterior margin of gh to middle of anal opening) and the total vaginal length (tvL).



Point	Description	Range of Values	Stage	Description
Aa (1)	Anterior vaginal wall 3cm proximal to the hymen	-3cm to +3cm	0	Aa, Ap, Ba, Bp = -3cm and C or D ≤ -(tvL -2)cm
Ba (2)	Most distal point of the remaining upper anterior vaginal wall	-3cm to +tvL		
C (3)	Most distal edge of cervix or vaginal cuff scar		1	Stage 0 criteria not met and leading edge < -1
D (4)	Posterior fornix (N/A if post-hysterectomy)			
Ap (5)	Posterior vaginal wall 3cm proximal to the hymen	-3cm to +3cm	2	Leading edge ≥ -1cm but ≤ +1cm
Bp (6)	Most distal position of the remaining upper posterior vaginal wall	-3cm to +tvL		
gh (7)	Hiatus genitalis, measured from middle of external urethral meatus to posterior midline hymen		3	Leading edge > +1cm but < +(tvL -2)cm
pb (8)	Perineal body, measured from posterior margin of gh to middle of anal opening			
tvL (9)	Total vaginal length (tvL), depth of vagina when point D or C is reduced to normal position		4	Leading edge ≥ +(tvL -2)cm

Endometrial Sampling Devices

Our endometrial samplers are designed for efficient, comfortable collection of endometrial tissue for histological evaluation, cancer screening, endometrial dating, or microbiological testing.

Each device features a smooth, flexible shaft and a precision-engineered tip to help clinicians obtain reliable specimens with minimal discomfort.

Preferred Curette

A family of single-use curettes available in multiple diameters and flexibilities to match different anatomical needs and clinical indications.



SKU	Description	Diameter	
BA8200	SureFlex Preferred Premium A premium 3.0 mm curette combining enhanced glide and controlled flexibility for improved patient comfort and superior clinician handling.	3.0mm	25 units/box
M0014	Standard Preferred Curette Designed for routine endometrial biopsy, the 3.0 mm diameter provides an ideal balance of patient comfort, tactile feedback, and consistent tissue sampling for everyday clinical use.	3.0mm	25 units/box
M0015	SureFlex Preferred Curette Also 3.0 mm in diameter, but with added flexibility to improve access in curved or challenging uterine anatomy while maintaining reliable sample collection.	3.0mm	25 units/box
M0016	EZ Preferred Curette With a larger 3.6 mm diameter, this model is designed for higher tissue yield when greater diagnostic volume is required.	3.6mm	25 units/box

Ensula™ Curette

A 3.0 mm suction curette with a Randall-type cutting edge, packaged with a twist-and-lock syringe. Designed for efficient collection of endometrial tissue for histological examination, especially useful when suction assistance is preferred.

SKU	Description	
RD003	Curette W/syringe	25 units/box



Cytology Brushes

Cytology brushes are used to collect samples from the cervix for PAP smears. The samples are examined to detect the early stages or the onset of cancer.


SKU	Description	Length	
TEA R1011	Cytology Brush	200mm	100units/pack



Dilators



Silicone Vaginal Dilator Sets



Each set includes 4 graduated sizes of dilators used to enlarge the vaginal for examination, surgical procedures, or relief of vaginismus. Single patient use only.

Large Set

Family	Length	OD Small	OD Large
DT-A/XL	150.0	30.0	38.0
DT-A/L	138.0	25.5	32.0
DT-A/M	128.0	22.5	28.0
DT-A/S	120.0	19.5	24.0

Medium Set

Family	Length	OD Small	OD Large
DT-B/XL	95.0	26.0	33.0
DT-B/L	86.0	22.5	28.0
DT-B/M	78.0	19.5	23.0
DT-B/S	75.0	12.5	17.0

Small Set

Family	Length	OD Small	OD Large
DT-C/XL	75.0	13.5	18.0
DT-C/L	65.0	12.0	16.0
DT-C/M	60.0	10.5	14.0
DT-C/S	50.0	9.5	13.0


Unit : mm

Bioteque Vaginal Dilator: CPT Code: 99070 Procedure Code: 57400

OS Finder Sets 25/Box

Bioteque offers disposable (blue) OS Finder set. The OS Finder set include an OS Finder, Canal Locator, and Fundus Dilator. The disposable OS Finders can be ordered separately. These sets are useful in locating and dilating the external and internal OS as well as the cervical canal.

SKU	Description	Length	Taper	O.D.
OS265BA-1	Disposable OSFinder	208mm	23mm	Ø 1.75 - 3.8mm
OS265BA-2	Disposable Canal Locator	208mm	45mm	Ø 1.75 - 3.8mm
OS265BA-3	Disposable Fundus Dilator	208mm	85mm	Ø 1.75 - 3.8mm



Uterine Sound 25/Box

Bioteque's latex-free Uterine Sounds are malleable, and disposable. Used for dilating and measuring the cervix. Graduated in centimeters with markings at 3cm, 5cm, 7cm, 9cm, and 11cm.

SKU	Description	Length	O.D.
BA100	Disposable Uterine Sound	248.5mm	3mm

Catheters

Word Bartholin Catheter Set

This latex-free Bartholin Catheter Set contains a silicone inflatable catheter, scalpel, syringe, and needle, for the treatment of Bartholin gland abscesses.

SKU	CPT Code	HCPT Code	
WC-400MM	56420	A4550	6units/pack

Uterine Manipulator Injector

Ideally designed for Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Minilaparotomy, Salpingoplasty, and Fertility Exams.

1. Pear shaped balloon for optimum hold and better seal at the internal OS
2. Easy insertion by reducing the size of cannula
3. Flexible sleeve keeps luer from cracking
4. No need for another instrument during the procedure

SKU	Description	O.D.	
UMI 4.5	Uterine Manipulator/ Injector	4.5mm	12units/box
UMI 4.5K	Uterine Manipulator/ Injector with clamp	4.5mm	12units/box

Uterine Injector

1. Sterile, disposable, ready to use
2. Curved Shaft (uterine sound shaped) for use in anteverted or retroverted uteri
3. Snap-on cervical stop custom set for every patient
4. Snap-on cervical stop V-shaped design for easy view of external cervical OS
5. Patient ambulatory once instrument properly inserted
6. Durable no-crack luer fitting on flexible collar for convenient injection
7. No spring or weights added

SKU	Description	Tube O.D.	
UI 4.0	Uterine Injector 4.0mm	4.0 mm	12units/box

HSG Stylet

SKU	Description	Length	
BA001	H/S Catheter Stylet	420mm	10units/box

HSG Catheter Sets Latex Free / with balloon

Bioteque’s latex-free HSG Catheter Sets, available in 5F &7F with a silicone balloon, are useful in evaluating the causes of abnormal bleeding, menstrual disorders, recurrent miscarriage, and unexplained infertility.

SKU	Description	
BA05F	5FrenchCatheters	10 units/ box
BA07F	7FrenchCatheters	10 units/ box
BA051F	5FrenchSureFlexCatheters	10 units/ box
BA071F	7FrenchSureFlexCatheters	10 units/ box

Bellota™ Sonohysterography Catheter

The balloonless Sono Catheter is a sterile disposable 5F catheter that consists of an insertion sheath and cervical cone. The adjustable cervical cone positioner helps with catheter placement while reducing leakage of saline at the external cervical OS.

SKU	Description	Length	
UIC5F	SonoCatheter w/o cervical cone 5Fr	327.5mm	25units/box
UIC5FA	SonoCatheter with cervical cone 5Fr	276.5mm	25units/box
UIC5FAX	SonoCatheter with cervical cone 5Fr	259.5mm	25units/box
UIC5FS	SonoCatheter/Uterine Introducing 5Fr	350.0mm	25units/box


CPT Code : 76831

Electrosurgical Devices

Bioteque offers a full range of electrosurgical devices designed for reliable performance and everyday clinical demands. Driven by quality, safety, and continuous improvement, dependable products and service support are provided for healthcare professionals worldwide.

Leep/Lletz Electrodes

Bioteque loop electrodes are designed for precise cervical tissue excision during LEEP/LLETZ procedures and are available in a wide range of shapes and sizes, including round, square, and specialized loop configurations to meet varying clinical requirements. The loop electrodes support accurate tissue removal while helping to limit unnecessary tissue excision. All electrodes are sterile, single-use, and individually sealed in procedure-ready blister packaging.



SKU	Description	
OB-1512	15mm x 12mm	Loop
OB-1515	15mm x 15mm	Loop
OB-2008	20mm x 8mm	Loop
OB-2010	20mm x 10mm	Loop
OB-2012	20mm x 12mm	Loop
OB-2015	20mm x 15mm	Loop
OB-2020	20mm x 20mm	Loop
OB-2510	25mm x 10mm	Loop
OB-3B	3mm	Ball
OB-5B	5mm	Ball

SKU	Description	
OB-1004S	10mm x 4mm	Square Loop
OB-0505	5mm x 5mm	Loop
OB-0505S	5mm x 5mm	Square Loop
OB-1005	10mm x 5mm	Loop
OB-1007	10mm x 7mm	Loop
OB-1008	10mm x 8mm	Loop
OB-1008S	10mm x 8mm	Square Loop
OB-1010	10mm x 10mm	Loop
OB-1010S	10mm x 10mm	Square Loop
OB-1507	15mm x 7mm	Loop
OB-1510	15mm x 10mm	Loop

Electrosurgical
ELECTRODE



Electrosurgical Pencil

SKU	Description	
OB-48020	Electrode Button Pencil	10units/pack

Electrosurgical Accessories


Bioteque offers a full range of LEEP/LLETZ accessories, including dispersive pads, tubing, and cables, designed to support safe and reliable electrosurgical procedures. All accessories are manufactured under strict quality standards to ensure consistent performance and compatibility in daily clinical use.



Smoke Tubing Set

6 Ft X 1/2 inch Smoke Evacuation Tubing
Smooth Bore Sets with Straw

SKU	Description	
OB-48030	Tubing Non-Sterile 6' long	10units/box



Smoke Tubing Set

Smoke Evacuation Tubing with straw
Fluted 6' Length, Closed Tube

SKU	Description	
OB-48034	Smoke Evacuation Tubing	10units/box



Ground pads

Ground pads with or without cable
Also fits Cooper Leep generator

SKU	Description	
OB-48050	Ground pads w/Cables	10units/pack
OB-48050C	Ground Pad w/Cables fits Cooper	10units/pack
OB-48055	Disposable Ground Pad w/Cable	5units/pack



ULPA Primary filter

Fits Stackhouse smoke evacuator

SKU	Description	
OB-48070	Prefilter for Smoke Evacuation Funnel	5units/pack
OB-48075	ULPA Primary Filter fits Stackhouse Minivac	10units/pack



Reducer

SKU	Description	
OB-48080	Reducer 3/8	10units/pack
OB-48082	Reducer 7/8	10units/pack



Monsel's Solution

Ferric Subulfate Solution Box/12 8ml bottles and 12 rayon tip applicators

SKU	Description	
9045055	Monsel Astrigen	12units/box

Instruments

Bioteque stainless steel instruments are manufactured from high-quality medical-grade stainless steel to ensure durability, corrosion resistance, and long-term performance in clinical environments. Designed for precision, balance, and reliable handling, each instrument supports consistent performance across a wide range of medical and surgical applications.

All instruments are suitable for routine cleaning, disinfection, and sterilization processes, ensuring safety, hygiene, and dependable reuse in daily clinical practice.

Biopsy Punches



Tischler-Morgan
SKU: BP018 10"



7.0 mm x 3.0 mm



Baby Tischler
SKU: BP019 10"
SKU: BP009 8.25"



4.2 mm x 2.3 mm



Burke
SKU: BP011 8.25"



5.5 mm x 3.25 mm



Tischler-Kevorkian
SKU: BP012 8.25"



9.5 mm x 3.0 mm



Kevorkian-Pacific
SKU: BP022 10"



8.0 mm x 3.0 mm



BP018



BP009



BP011


*bite shape not to scale, for reference only

SKU	Description
BP008	Tischler PCH Sur-Lok 8.25
BP009	Baby Tischler Bitesize 4.2 X 4.3mm
BP011	Oval PCH Sur-Lok 8.25
BP012	Tisch-Kev Punch Sur-Lok 8.25
BP013	Biopsy Punch Kevorkian (pacific) 8.25
BP017	Sullivan PCH 8.25


SKU	Description
BP018	Tischler PCH Sur-Lok 10
BP019	Baby Tischler PCH Sur-Lok 10
BP020	Oval Biopsy Punch 10
BP021	Tischler-Kevorkian Sur-Lok Pch 10
BP022	Kevorkian-Pacific PCH Sur-Lok 10
BP023	Baggish Sur-Lok Pch 10



BP111



BP103



BP105

	Bite Size	Biopsy Punch	8.25" SKU	10" SKU
	4.2 mm x 2.3 mm	Roto Baby Tischler	BP111	BP104
	7.0 mm x 3.0 mm	Roto Tischler	BP109	BP103
	5.5 mm x 3.25 mm	Roto Oval (Burke)	N/A	BP105
	7.0 mm x 3.0 mm	Roto Tischler Morgan	BP114	N/A

*bite shape not to scale, for reference only



Foerster Forcep Straight 9.5" (SKU: FC008)



Bozeman Uterine Forcep CVD 10 " (SKU:FC206)

Forceps

SKU	Description
FC006	I.U.D. Removal FCP 12"
FC008	Foerster Forcep Straight 9.5 "
FC106	Schroder Tenaculum Forceps 9.5"
FC107	Braun Tenaculum Forceps 9.5"

Speculums Series

Lletz Speculums & Instruments (Autoclavable)



Lletz Graves Vaginal Speculum
w/ Evacuation Tube
Large
SKU: LZ003



Lateral Retractor
4.8cm

SKU: LZ202



Lletz Thumb Dress Forcep
8" SKU: LZ101



Lletz Tenaculum Forcep 9.5"
SKU: LZ108

SKU	Description
LZ001	Lletz Graves Speculum SM E/T 3"
LZ002	Lletz Graves Speculum Med E/T w/EVAC
LZ003	Lletz Graves Speculum Large E/T w/EVAC
LZ004	Lletz Graves Speculum E/T XL 7"
LZ006	Lletz Graves-Snowman Speculum E/T 6"
LZ007	Lletz Pederson Speculum SM E/T 3"
LZ008	Lletz Pederson Speculum MD E/T 4"
LZ009	Lletz Pederson Speculum LG E/T 4"
LZ010	Lletz Pederson Speculum XL E/T 7"

SKU	Description
LZ011	Lletz Speculum Weis Graves Med RT W/EVAC
LZ012	Lletz Speculum Weis Graves Large RT W/EVAC
LZ023	Lletz W-PEDS R XL E/T 7"
LZ024	Lletz View-More Speculum E/T SM 3"
LZ025	Lletz View-More Speculum E/T MD 4"
LZ026	Lletz View-More Speculum E/T LG 4.75"
LZ028	Lletz View-Maxi Speculum E/T GR LG 4.75"
LZ029	Lletz View-Maxi Speculum E/T MD 4"

SKU	Description
LZ034	Lletz Ultra-view MD E/T 4"
LZ041	Lletz Kogan Speculum 4-8mm 11
LZ101	Lletz Thumb dress Forcep 8"
LZ108	Lletz Tenaculum Forcep 9.5"
LZ202	Lateral Retractor 4.8cm
LZ203	Lletz Lateral Retractor V-M 10
LZ205	Lletz R.W. Retractor Med BL 7.75
LZ206	Lletz R.W. Retractor LG BL 7.75




Standard 30mm opening



View-More 40mm opening



View-Maxi 70mm opening



Graves Small
Standard 30mm Opening
SKU: SP001

Moore Graves
Standard 30mm Opening
SKU: SP005

Graves Medium
Standard 30mm Opening
SKU: SP002

Graves Large
Standard 30mm Opening
SKU: SP003
View-More 40mm Opening
SKU: SP216
View-Maxi 70mm Opening
SKU: SP218



Weisman-Graves Vaginal Speculum
Medium
R/Open
SKU: SP201



Pederson Vaginal
Speculum Large
SKU: SP105



Graves Vaginal Speculum
Large
SKU: SP003

SKU	Description
SP001	Graves Vaginal Speculum SM 3
SP002	Graves Vaginal Speculum MD 4
SP003	Graves Vaginal Speculum LG 4.5
SP004	Graves Vaginal Speculum XXL 7"
SP005	Moore Graves Speculum 4
SP006	Graves-Snowman Speculum 6

SKU	Description
SP101	Pederson VAG Speculum XS 2.25
SP102	Pederson VAG Speculum EX NAR 4
SP103	Pederson VAG Speculum SM 3
SP104	Pederson VAG Speculum Medium
SP105	Pederson VAG Speculum LG 4.5
SP106	Pederson VAG Speculum XXL 7"

SKU	Description
SP201	Weisman Speculum Graves R Med 4
SP202	Weisman Speculum Graves R LG 4.5
SP215	View-More Speculum Med 4
SP217	View-Maxi Graves Speculum Med 4
SP218	View-Maxi Graves Speculum LG 4.5
SP219	View-Maxi PEDS Speculum Med 4

Stainless Steel Curette

Novak Endocervical Curette



Novak



Kevorkian

SKU	Description
CT001	Randall Endometrium Biopsy Suction Curette with Luer-Lock 4mm 9.75"
CT002	Novak Endometrium Biopsy Suction Curette with Luer-Lock 1mm 9.75"
CT003	Novak Endometrium Biopsy Suction Curette with Luer-Lock 2mm 9.75"
CT004	Novak Endometrium Biopsy Suction Curette with Luer-Lock 3mm 9.75"
CT005	Novak Endometrium Biopsy Suction Curette with Luer-Lock 4mm 9.75"
CT006	Kevorkian Curette without Basket 12
CT007	Kevorkian Curette with Basket 12
CT008	Townsen Curette without Basket 12



Sims Uter CUR#00 CT101

Sims Uter CUR#0 CT102

Sims Uter CUR#1 CT103

Sims Uter CUR#2 CT104

Sims Uter CUR#3 CT105



3-4mm
5-6mm
7-8mm
9-10mm
11-12mm
13-14mm
15-16mm
17-18mm

Hegar Dilator



9/10
11/12
13/14
15/16
17/18
19/20

Hank Dilator



13/15
17/19
21/23
25/29
33/35
37/39
41/43

Pratt Dilator

Stainless Steel Dilators

SKU	Description
DL003	Mini Dilator Set 1.0, 1.5, 2.0, 2.5, and 3.0mm
DL101	Dilator Cervical Set of 6, 1-6 mm
DL201	Hegar Dilator Set of 8, single / double ended
DL300	Hank Dilator Double Ended 7 / 8 Fr
DL301	Hank Dilator Set of 6, 11"
DL302	Hank Dilator 9 / 10 Fr, 11"
DL303	Hank Dilator 11 / 12 Fr, 11"

SKU	Description
DL304	Hank Dilator 13 / 14 Fr, 11"
DL305	Hank Dilator 15 / 16 Fr, 11"
DL306	Hank Dilator 17 / 18 Fr, 11"
DL307	Hank Dilator 19 / 20 Fr, 11"
DL401	Pratt Uterine Dilator Set of 8, double ended 11.5"
DL501	Denniston Dilator Set of 5

Pessary Selection Chart

Pessaries may be used to relieve the following symptoms

Type/ symptoms		 Prolapse (1-2nd Degree)	 Prolapse (2-3 Degree)	 Cystocele	 Rectocele	 Stress Incontinence	 Retro Displacement
	Ring	✓		✓			
	Cube		✓	✓	✓		
	Gellhorn		✓	✓			
	Short-stem Gellhorn		✓	✓			
	Donut		✓	✓	✓		
	Shatz	✓	✓				
	Gehring		✓	✓	✓		
	Dish	✓	✓	✓		✓	
	Cup	✓	✓	✓		✓	
	Oval	✓		✓			
	Ring with Knob	✓		✓		✓	
	Marland	✓	✓	✓		✓	
	Hodge	✓		✓		✓	✓
	Incontinence Ring					✓	
	FlexiShelf		✓				

RETURN POLICY

Product(s) may be returned at Bioteque's expense for the following reasons only:

- The received product(s) is defective.
- The product(s) is incorrect.
- The product(s) was not ordered.

Bioteque will credit and accept returns of non-sterile product(s) at a 25% restocking fee only if:

- The product(s) is returned in the ORIGINAL packaging and marketable condition and
- A written reason for the return is provided and
- A copy of Bioteque's packing list or invoice is included

A Return Merchandise Authorization (RMA) number must be obtained from Bioteque by calling (408) 526- 9398.

The RMA number must appear on the outside of the package.

All returns must be sent to Bioteque America, Inc.

1302 Exchange Drive Suite #150 Richardson, TX 75081

NO RETURNS WILL BE ACCEPTED AFTER 90 DAYS OF RECEIPT. THERE ARE NO REFUNDS OR CREDIT FOR THE RETURN OF STERILE PRODUCTS.

FOR ANY QUESTIONS REGARDING THE RETURN POLICY, CONTACT BIOTEQUE AT 408-526-9398

ORDERING INFORMATION

Distribution Center:

1302 Exchange Drive Suite #150 Richardson, TX 75081

Remittance To:

Bioteque America Inc.

1302 Exchange Drive Suite #150 Richardson, TX 75081

Phone Orders To:

+1 (800) 889-9008

+1 (408) 526-9398

Sales Inquiries To:

sales@bioteque.com


General information/Questions and Concerns:

Info@bioteque.com

customerservice@bioteque.com

Credit Card Orders:

Visa, Mastercard, AMEX



With 30-years experience and global operations, we serve hospitals, clinics, and distributors worldwide.

Bioteque America Inc.
1302 Exchange Drive
Suite 150
Richardson, TX 75081

Bioteque Surgical Europe Ltd.
Unit 1E, Three Rock Road,
Sandyford Business Park
Dublin 18, Ireland

Bioteque Taiwan Ltd.
13F., No. 9, Ln. 130, Sec. 3,
Minsheng E. Rd.
Songshan Dist., Taipei City

Bioteque Philippines
2nd Flr. Unit 1 GTSRR Commercial Bldg.
Bahak Pobiacion, 6002
Liloan, Cebu Philippines

Product Catalog
V 6.0
Dec 2025

