

Pessaries

Bioteque offers a complete range of medical-grade silicone pessaries in a variety of shapes and sizes to accommodate different patient needs. Each product is manufactured under strict quality control systems and complies with FDA, cGMP, CE, QMS, MDD, ARTG, and ISO 13485 standards, ensuring safety, consistency, and reliability.



Ring

The Ring pessary, available with or without a supporting membrane, is primarily used for mild first- degree uterine prolapse and cystocele.

SKU (W/O support)	R#0	R#1	R#2	R#3	R#4	R#5	R#6	R#7	R#8	R#9
Dimensions (mm)	44	50	57	63	69	75	82	88	95	101
Dimensions (inch)	1.75	1.95	2.25	2.50	2.70	2.95	3.25	3.45	3.75	4.00
SKU (with support)	Rs#0	Rs#1	Rs#2	Rs#3	Rs#4	Rs#5	Rs#6	Rs#7	Rs#8	Rs#9



Ring with Knob

The Ring with Knob pessary, available with or without a supporting membrane, is primarily used for mild first-degree uterine prolapse, cystocele, and stress urinary incontinence.

SKU (W/O support)	RK#0	RK#1	RK#2	RK#3	RK#4	RK#5	RK#6	RK#7
Dimensions (mm)	44	50	57	62.5	69	75	82	87.5
Dimensions (inch)	1.75	1.95	2.25	2.45	2.70	2.95	3.25	3.45
SKU (with support)	RKs#0	RKs#1	RKs#2	RKs#3	RKs#4	RKs#5	RKs#6	RKs#7



Incontinence Ring

The Incontinence Ring is primarily used for stress urinary incontinence, mild first-degree uterine prolapse, and cystocele.

SKU	IRO	IR1	IR2	IR3	IR4	IR5	IR6	IR7	IR8	IR9
Dimensions (mm)	44	51	57	64	70	76	83	89	95	102
Dimensions (inch)	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00



*Bioteque Exclusive

Oval

The Oval pessary, available with or without a supporting membrane, is primarily used for mild first-degree uterine prolapse and cystocele. The pessary's shape helps to limit rotation.

SKU(W/O support)	OV#1	OV#2	OV#3	OV#4	OV#5	OV#6	OV#7	OV#8	OV#9
Dimensions (mm)	50.5	57	63.5	69.5	76	83	89	96	101
Dimensions (inch)	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.80	4.00
SKU (with support)	OVs#1	OVs#2	OVs#3	OVs#4	OVs#5	OVs#6	OVs#7	OVs#8	OVs#9



*Bioteque Exclusive

Marland

The Marland pessary, available with or without a supporting membrane, is primarily used for second or third-degree uterine prolapse, cystocele, rectocele. The pessary's elevated platform may also relieve stress urinary incontinence.

SKU (W/O support)	M#2	M#3	M#4	M#5	M#6	M#7	M#8
Dimensions (mm)	58	64	70	76	82	88	94
Dimensions (inch)	2.30	2.50	2.75	3.00	3.25	3.45	3.70
SKU (with support)	Ms#2	Ms#3	Ms#4	Ms#5	Ms#6	Ms#7	Ms#8



Gehrung

The malleable Gehrung pessary with a supporting membrane is primarily used for significant cystocele or rectocele, and second or third-degree uterine prolapse. The pessary can be manually shaped to accommodate individual anatomies. Use caution and follow doctor's instructions before undergoing MRI.

SKU	GH#0	GH#1	GH#2	GH#3	GH#4	GH#5	GH#6	GH#7	GH#8	GH#9
Dimensions (mm)	50	55	60	65	70	75	80	85	90	95
Dimensions (inch)	1.95	2.15	2.35	2.55	2.75	2.95	3.15	3.35	3.55	3.75



Shaatz

The Shaatz pessary is primarily used for second or third-degree uterine prolapse and cystocele. Its convex shape ensures a snug fit drainage holes aid in removal.

SKU	SH#0	SH#1	SH#2	SH#3	SH#4	SH#5	SH#6	SH#7	SH#8	SH#9
Dimensions (mm)	38	45	51.5	57	63	70	76	82	88	95
Dimensions (inch)	1.50	1.75	2.05	2.25	2.50	2.75	3.00	3.25	3.45	3.75



Donut

The Donut pessary is primarily used for second or third-degree uterine prolapse, cystocele, and rectocele. This pessary, often used by post-menopausal women, is easily compressed for insertion and removal.

SKU	D#0	D#1	D#2	D#3	D#4	D#5	D#6	D#7
Dimensions (mm)	50	59	65	69	76	82	88	95
Dimensions (inch)	1.95	2.30	2.55	2.70	3.00	3.25	3.45	3.75

Bioteque Pessary CPT Code: A4562 Procedure Code: 57160



*Bioteque Exclusive

FlexiShelf

The Flexishelf pessary is primarily used for second or third-degree uterine prolapse. It has a unique thin curved-base design, rigid inner layer and soft outer layer. The curved cone-shaped handle aids in removal.

SKU	FS06	FS08	FS10	FS12
Dimensions (mm)	57	64	70	76
Dimensions (inch)	2.25	2.50	2.75	3.00



Cube

The Cube pessary, available with or without drains, is primarily used for cystocele, rectocele, second or third-degree uterine prolapse, and complete procidentia.

SKU (w/O drains)	CU#0	CU#1	CU#2	CU#3	CU#4	CU#5	CU#6	CU#7	CU#8	CU#9	CU#10
Dimensions (mm)	25	29	33	37	41	45	50	56	63	70	75
Dimensions (inch)	1.00	1.15	1.30	1.45	1.60	1.75	1.95	2.20	2.50	2.75	2.95
SKU (with drains)	CUD#0	CUD#1	CUD#2	CUD#3	CUD#4	CUD#5	CUD#6	CUD#7	CUD#8	CUD#9	CUD#10



Gellhorn

The Gellhorn pessary with drains is primarily used for cystocele, rectocele, and second or third-degree uterine prolapse. This pessary's stem with knob is easily folded to aid with insertion and removal.

SKU (with drains)	Gd#0	Gd#1	Gd#2	Gd#3	Gd#4	Gd#5	Gd#6	Gd#7	Gd#8	Gd#9	Gd#10
Dimensions (mm)	38	45	51.5	57	63	70	76	82	88	95	101.6
Dimensions (inch)	1.50	1.75	2.05	2.25	2.50	2.75	3.00	3.25	3.45	3.75	4.00



Short Stem Gellhorn

The short-stem Gellhorn pessary with drains is primarily used for cystocele, rectocele, and second or third-degree uterine prolapse. It provides a better fit for some women,

SKU (with drains)	GS#0	GS#1	GS#2	GS#3	GS#4	GS#5	GS#6	GS#7	GS#8	GS#9
Dimensions (mm)	38	44	51	57	63.5	70	76	83	89	95
Dimensions (inch)	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75



Dish

The Dish pessary, available with or without a supporting membrane, is primarily used for second or third-degree uterine prolapse, cystocele, rectocele and stress urinary incontinence.

SKU (W/O support)	DSH#0	DSH#1	DSH#2	DSH#3	DSH#4	DSH#5	DSH#6	DSH#7
Dimensions (mm)	48.5	53.5	59	64	69	75	79	84
Dimensions (inch)	1.90	2.10	2.30	2.50	2.70	2.95	3.10	3.30
SKU (with support)	DSHs#0	DSHs#1	DSHs#2	DSHs#3	DSHs#4	DSHs#5	DSHs#6	DSHs#7



Cup

The Cup pessary, available with or without a supporting membrane, is primarily used for first or mild second-degree uterine prolapse. This pessary is similar to the dish without a knob.

SKU (W/O support)	CP#0	CP#1	CP#2	CP#3	CP#4	CP#5	CP#6	CP#7
Dimensions (mm)	48.5	53.5	59	64	69	74	79	84
Dimensions (inch)	1.90	2.10	2.30	2.50	2.70	2.90	3.10	3.30
SKU (with support)	CPs#0	CPs#1	CPs#2	CPs#3	CPs#4	CPs#5	CPs#6	CPs#7



Hodge

The malleable Hodge pessary, available with or without a supporting membrane, is primarily used for first or mild second-degree prolapse and cystocele. This pessary, is usually prescribed for reproductive-aged women who prefer conservative therapy. Use caution and follow doctor's instructions before undergoing MRI.

SKU (W/O support)	HD#0	HD#1	HD#2	HD#3	HD#4	HD#5	HD#6	HD#7	HD#8	HD#9
Dimensions (mm)	65	70	75	80	85	90	95	100	105	110
Dimensions (inch)	2.55	2.75	2.95	3.15	3.35	3.55	3.75	3.95	4.15	4.35
SKU (with support)	HDs#0	HDs#1	HDs#2	HDs#3	HDs#4	HDs#5	HDs#6	HDs#7	HDs#8	HDs#9

Silicone Inflatable Donut



*Bioteque Exclusive

The Inflatable donut is primarily used for mild second or third degree cystocele in combination with procidentia/prolapse.

SKU	Dimensions	Dimensions
ID 2	51mm	2.00in
ID 3	57mm	2.25in
ID 4	64mm	2.50in
ID 5	70mm	2.75in