



PRODUCT LITERATURE | PESSARIES

Pessary Selection Chart

Official Manufacturer Selection chart | English

The following pages reproduce the manufacturer's official document without alteration. For availability, quotations and support in Canada, please contact Beauray Inc.

CANADIAN SALES & SUPPORT

Beauray Inc.

info@beauray.co

+1 437 607 5116 | Thornhill, Ontario, Canada

For Canadian product inquiries, contact Beauray Inc.

Pessary Selection Chart

Pessaries may be used to relieve the following symptoms

Type/ symptoms		 Prolapse (1-2nd Degree)	 Prolapse (2-3 Degree)	 Cystocele	 Rectocele	 Stress Incontinence	 Retro Displacement
	Ring	X		X			
	Cube		X	X	X		
	Gellhorn		X	X			
	Short-stem Gellhorn		X	X			
	Donut		X	X	X		
	Shaatz	X	X				
	Gehrung		X	X	X		
	Dish	X	X	X		X	
	Cup	X	X	X		X	
	Oval	X		X			
	Ring with Knob	X		X		X	
	Marland	X	X	X		X	
	Hodge	X		X		X	X
	Incontinence Ring					X	
	FlexiShelf		X				